

# editsetting.ini

Setting		Description
<b>[Version]</b>		<b>CDF Header</b>
sLab = bouma		Data lab origin
sVersion = bouma_unknown		Version of data (metadata namely, .xml, .dat, .ini, etc.)
sState = struct+angio+ps		Processing state (N/A 2021/01/07) *See End
<b>[CplxTomSetting]</b>		<b>CDF Header</b>
nSamples = 1688		Number of samples acquired per a-line
nLines = 2048		Number of lines in the struct dataset
nFrames = 28		Number of b-scans/frames
nZscans = 2048		Number of z points in the output
nZoomLevel = 2		Post-processing zoom level
nZoomLevelRT = 2		Real-time zoom level
nIndexLow = 0		Cropping index low (for projections and hsv output)
nIndexHigh = 0		Cropping index high (for projections and hsv output)
nFrameInterval = 1		Process every X frames, this interval
bFlipaline = 1		Flip Aline at after tomogram reconstruction
sMappingfilename =		Mapping filename (N/A 2021/01/07)
sDispersionfilename =		Dispersion filename (N/A 2021/01/07)
sConfigfilename =		Config filename (N/A 2021/01/07)
sScanpatternfilename =		Scan pattern filename (N/A 2021/01/07)
afDemodulation = 0.5,0,1,0,0,0,		Demodulation ratio
fClockRateMHz = 85		Clock rate of the acquisition card
<b>[StructTomSetting]</b>		<b>CDF Header</b>
fReflow = -40		Lower limit of log struct frame before uint8
fRefhigh = 130		Upper limit of log struct frame before uint8
bInvertgray = 0		Invert angio greyscale
<b>[AngioTomSetting]</b>		<b>CDF Header</b>
fReflow = -40		Lower limit of log weight frame before uint8
fRefhigh = 130		Upper limit of log weight frame before uint8
bInvertgray = 0		Invert angio greyscale
<b>[PsTomSetting]</b>		<b>CDF Header</b>
nMaxRet = 100		Maximum retardance, used to normalization over image volumes
nBinFract = 3		Fraction of spectrum for spectral bin (3=5bins,4=9bins etc.)
nZOffset = 10		Delta Z for local birefringence
fZResolution = 5		Resolution of system in Z

## editsetting.ini

nXFilter = 11		Across Aline (X) filter
nZFilter = 1		Along Aline (Z) filter
nOpFilter = 1		Out of plane frames to average (if OOP option is True)
fDopThresh = 0.7		dop threshold for PS contrast
nThetaOffset = 0		Theta offset (N/A 2021/01/07)
<b>[HSVSetting]</b>		<b>CDF Header</b>
nThetaRef = 0		Theta offset (rotates colormap)
bHueCCW = 0		(N/A 2021/01/07)
fOpacity = 0.01		Opacity of HSV projection
nHsvCropLow = 0		(N/A 2021/01/07) -> Use [CplxTomSetting] nIndexLow/
nHsvCropHigh = 0		(N/A 2021/01/07) -> Use [CplxTomSetting] nIndexHigh
nDopWeightLow = 20		Contrast weight from dop
nDopWeightHigh = 130		Contrast weight from dop
nStructWeightLow = 30		Contrast weight from struct
nStructWeightHigh = 100		Contrast weight from struct
nRetWeightLow = 10		Contrast weight from ret
nRetWeightHigh = 100		Contrast weight from ret
nMaskThresholdsDOP = 30		Lower threshold for dop Mask
nMaskThresholdsRet = 30		Lower threshold for ret Mask
nMaskThresholdsStruct = 30		Lower threshold for struct Mask
<b>[ProcOptions]</b>		<b>CDF Header</b>
bOOPAAveraging = 0		Perform out of plane averaging
bFastProcessing = 0		Perform fast downsampling for bi-seg scans
bSpectralBinning = 1		Perform spectral binning in PS reconstruction
bCorrectSystemOA = 0		Correct the system's optic axis, create "symmetric" conditions
bCorrectSystemDiat = 0		Correct system diattenuation (N/A 2021/01/07)
nNFramesOACorr = 10		Number of equally spaced frames within volume to do correction **Input can be integer or a list: ie. 110,121,133,150
bComputeBackground = 0		Compute median based background from dataset
nNFramesBGCorr = 10		Number of equally spaced frames within volume to do correction **Input can be integer or a list: ie. 110,121,133,150
bRotCartesianOutput = 0		Output circular image volumes for endoscopy (N/A 2021/01/07)
bMaskOutput = 0		Mask ps output using dop/struct/ret
bGenerateProjections = 0		Generate projections
sProjState = 'struct+angio'		Frame types to be projected
sProjType = 'max+mean+std+sum'		Type of projection to be performed

\*sState options include: 'tomo+kspace+stokes+struct+angio+ps+**hsv+oa**+'  
       Signifies datatypes that are *processed but not written*, including these states

in the string will cause the processor to write them to the Processed directory.